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
English

INTENDED USE

Medical device used to connect a breathing circuit or resuscitator to an Endotrachealtube or Face Mask, to enhance mobility at the patient end.

INDICATIONS

- To enhance mobility of the breathing circuit at the patient end during anaesthesia or ventilation.

CONTRAINDICATIONS

- Relative contraindication: ventilation of children and infants.

Further contraindications are not known.

SAFETY INFORMATION

- Read and follow the instructions for use carefully before using the product.
- The product must only be used by trained medical personnel.
- When using connecting tubing with a length of 30 or 45 cm, the silicone ring must be removed before application.
- Ensure correct and tight connection.
- For disconnection of stuck connectors, the manufacturer offers an appropriate uncoupling wedge.
- The product is suitable for conveyance of oxygen, air, anaesthetic gases and their gas mixtures. Anaesthetic gases used must not interact with the product materials (see chapter "Material Data").
- Due to a possibility to cause static charge, the product must not be used with or in striking distance of flammable anaesthetic gases.
- The user must pay attention for efficient ventilation during application.
- Prior to use, the products must be inspected visually for damages (cracks, breakage, etc.) and a functional check must be performed (see chapter "Functional check").
- The product is for single use and must not be reprocessed. Reprocessing influences the function of the product. Any reuse involves the potential risk for cross contamination.
- The product must not be used if the packaging is damaged or the shelf life has elapsed.

TECHNICAL DATA

Patient Group	Tube length	Intended delivered volume	Compliance (ml/cmH ₂ O)	Flow resistance (cmH ₂ O/l/min)
Adult	15 cm	≥ 300 ml	C@60hPa: 0,20 ml/hPa	R@30l/min: 0,0085 hPa/l/min
	30 cm	≥ 300 ml	C@60hPa: 0,30 ml/hPa	R@30l/min: 0,0092 hPa/l/min
	45 cm	≥ 300 ml	C@60hPa: 0,28 ml/hPa	R@30l/min: 0,0103 hPa/l/min

FUNCTION CONTROL

Prior to use check the products for foreign bodies and patency. All connectors must be connected to the tubing.

Check the system for leaks as follows:

Connect the Catheter Mount Tubing to the breathing system and seal the patient end. Use a pressure of 60 cmH₂O and check if the pressure drops or if air escapes within one minute.

Damaged products must not be used.

MATERIAL DATA

Tube	Polyethylene (PE)
Connector	Polycarbonate (PC)
Bronchoscope / Suction cap	Silicone

SHELF LIFE

Expiry date: see product label.

STORAGE AND TRANSPORT CONDITIONS

- Protect from heat and store in a dry place.
- Keep away from sunlight and light sources.
- Store and transport in the original packaging.

DISPOSAL

Used or damaged products must be disposed acc. to applicable national and international legal directives.

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