

RTR II

Rectal Tampon by Roche

Rectal Tampon by Roche


Inflatable tampon to stop bleeding, for pain control and to release gases and fluids after anal surgery for adults.

Patients after anal surgery are affected by pain and bleeding. The classic anal medications are generally uncomfortable, glue with the blood and are difficult and painful to remove.

The RTR II is more comfortable for the patient, stops bleeding, allows post-operative pain control and pain free removal of the device. The inner tube allows the release of gases and fluids.

Indications: Surgical procedures in the anal canal

- Haemorrhoids
- Anal fistula
- Anal fissures
- Condyloma
- Anal tumours
- Anal polyps



Order information

**Rectal Tampon RTR II by Roche
with inner tube**

sterile, for single use

Length: 85 mm, Material: Silicone


REF 39-20-200-1

Box 10

VBM Medizintechnik GmbH

Einsteinstrasse 1 | 72172 Sulz a. N. | Germany

Tel.: +49 7454 9596-0 | Fax: +49 7454 9596-33 | e-mail: info@vbm-medical.de | www.vbm-medical.com



0123

P317-2.0GB/05.16