


Stylets, Introducers & Tube Exchangers

S-Guide

Malleable intubating stylet for video laryngoscopy

- for single use, sterile
- lumen-design, O₂ Connector included
- graduation marks
- unique design:
 - 1 soft tip
 - 2 malleable segment (up to length 42 cm)
 - 3 flexible segment (from length 42 cm)

Size	Length	for ET Tube	REF	Box
15 Fr	65 cm	≥ 6.0 mm	33-90-650-1	5



METTS (Muallem ET Tube Stylet)

Malleable intubation stylet

- core made of metal, malleable (maintains curvature)
- flexible preformed tip
- graduation marks
- for single use, sterile

Size	Length	for ET Tube	REF	Box
8 Fr	40 cm	≥ 3.5 mm	33-08-400-1	5
12 Fr	65 cm	≥ 5.0 mm	33-12-650-1	5
14 Fr	65 cm	≥ 6.0 mm	33-14-650-1	5



Introducer

To facilitate endotracheal intubation

- lumen-design, O₂ Connector included
- rigid preformed tip
- graduation marks
- for single use, sterile

Size	Length	for ET Tube	REF	Box
15 Fr	65 cm	≥ 6.0 mm	33-70-650-1	5




METTI (Muallem ET Tube Introducer)

To facilitate endotracheal intubation

- core made of plastic, semi-rigid (flexible)
- flexible, preformed tip
- graduation marks
- for single use, sterile

Size	Length	for ET Tube	REF	Box
12 Fr	80 cm	≥ 5.0 mm	33-12-800-1	5
14 Fr	80 cm	≥ 6.0 mm	33-14-800-1	5



Pocket Introducer

Ideal space solution for emergency bags and vehicles, helicopters, lab coat pockets, crash carts, etc.

- folded to 20 cm, unfolds to 65 cm
- ready for use within seconds
- preformed tip
- graduation marks
- for single use, sterile

Size	Length	for ET Tube	REF	Box
15 Fr	65 cm	≥ 6.0 mm	33-80-650-1	5

in robust plastic cartridge

Size	Length	for ET Tube	REF	Box
15 Fr	65 cm	≥ 6.0 mm	33-81-650-1	1



Tube Exchanger

For extubation and exchange of ET Tubes

- lumen-design, O₂ Connector included
- graduation marks
- for single use, sterile

Size	Length	for ET Tube	REF	Box
11 Fr	80 cm	≥ 4.0 mm	33-45-800-1	5
14 Fr	80 cm	≥ 5.0 mm	33-46-800-1	5
19 Fr	80 cm	≥ 7.0 mm	33-47-800-1	5



Stylet

To preform the ET Tube

- soft, atraumatic tip
- with and without moveable silicone connector
- reusable

Size	Length	for ET Tube	without silicone connector	Box	with movable silicone connector (15 mm I.D.)	Box
small	34 cm	2.5 - 3.5 mm	REF 90-10-111	25	REF 90-10-001	10
medium	38 cm	4.0 - 4.5 mm	REF 90-10-222	25	REF 90-10-002	10
large	45 cm	5.0 - 6.5 mm	REF 90-10-333	25	REF 90-10-003	10
large	45 cm	7.0 - 11.0 mm	REF 90-10-444	25	REF 90-10-004	10

