

Sulz am Neckar, November 2015

The new Surgicric

3 techniques – 3 cricothyrotomy sets for securing airways

Lancing of the airway in the vicinity of the larynx is a life-saving technique in the context of difficult airway management. This surgical airway is the ultima ratio measure of a „cannot intubate – cannot ventilate“ scenario according to all airway algorithm – in preclinical cases, in emergency situations, in intensive care as well as in surgery.

Various techniques for this life-saving emergency measure exist such as the surgical cricothyrotomy and the percutaneous cricothyrotomy.

Surgicric is a cricothyrotomy set for maintaining ventilation during an obstruction of the upper airways. VBM offers three pre-assembled set variations containing the necessary components for the conduct of the three main cricothyrotomy techniques. The heart of these sets is the special combination of tracheal tubes and dilatator. The locking mechanism and the soft point maximize patient safety and lower the risk of injury. The height-adjustable flange enables the individual adjustment of the tube to the anatomy of the patient.

VBM provides a cricothyrotomy set for the conduct of the Rapid Four-Step Technique by offerings Surgicric I. Surgicric II is applied according to the classic surgical technique. Surgicric III enables a cricothyrotomy according to the Seldinger Technique.

A further advantage of the sets is the unique instrument bag. It offers the user a clear overview of all components as well as enabling sterile use in any surrounding. The small pack size of the instrument bag is ideal for use in emergency bags.

Further information is available on the website: www.vbm-medical.com

About VBM Medizintechnik GmbH:

VBM Medizintechnik is headquartered in Sulz am Neckar (Germany) and employed approximately 200 persons. The private owned company develops and produces innovative device for Airway Management, Accessories for Anaesthesia and intensive Care as well as tourniquet systems for surgical operations in a bloodless field. VBM devices are present in more than 100 countries throughout distribution partners or own subsidiaries.

Contact:

Phone +49 7454 9596 0
presse@vbm-medical.de

Related keywords

Airway securing, Airway Management, difficult intubation, endotracheal Intubation

**Image**

File name: **VBM_Surgicric_Tubus_Dilatator.jpg**

Printing of these images is free of charge.
Please quote „VBM Medizintechnik“ as source.